





10611007
Aether 8 Eco
Vacuum Pump
User Manual



lable	e of Contents		
I. Disc	claimer & Warranty	3	
LIMITE	ED WARRANTY & DISCLAIMER	3	
LIMITA	ATION OF LIABILITY	6	
II. Safe	ety	7	
2.1	After Sales Support	9	
III. Aet	ther 8 Eco Vacuum Pump	10	
3.1 Co	omponents	10	
3.2	Before using the vacuum pump	11	
3.3	Post Processing Vacuum Pump Maintenance	11	
IV. Ted	chnical Specifications	12	
4.1	Specifications	12	
4.2	Expanded Drawing	13	
VI. Pro	oper Use	14	
VII. Ins	spection	14	
The package includes the following items			
VIII. M	laintenance	15	
8.1	Vacuum pump oil change procedure	15	
IX. Tro	oubleshooting	16	
9.1	Failure to start	16	
9.2	Oil leakage	16	
9.3	Failure to attain vacuum	16	



### I. Disclaimer & Warranty

- The intellectual property rights for this manual belongs to Holland Green Science Corporation (hereinafter referred to as "the Company").
- Products of the company are under the patent protection of China and other countries and regions (inclusive of obtained patents and patents currently under application).
- The company reserves the right to change product specifications without notice.
- The contents of the manual are subject to change without notice.-The manual's content
  must not be copied, extracted, or modified in any form without the Company's
  permission.

#### LIMITED WARRANTY & DISCLAIMER

Products sold by the Company are intended solely for the use(s) specified by the manufacturer and are not to be used for any other purposes including, but not limited to, unauthorized commercial purposes. Buyer warrants and represents that Buyer will properly test, use, qualify, and/or validate each Product for intended use, and, to the extent authorized, manufacture and market any final articles made from the Company's Products in accordance with the practices of a reasonable person who is an expert in the field, and in strict compliance with all applicable national, state and local laws and regulations. Buyer acknowledges and understands that the purchase of Products only conveys to Buyer the non-transferable right for only the Buyer to use the Products purchased in compliance with the applicable manufacturer's Published User Statement, Limited User Statement, or Limited License, if any. While this manual contains certain representations of the potential dangers from the instrument or equipment, Buyer is solely responsible for conducting any research necessary to learn the hazards involved for any of its intended uses of the instrument or equipment purchased from the Company and to properly warn its customers, employees and other personnel who may be exposed to such instrument or equipment of any risks involved in using or handling the instrument or equipment.

The Company warrants that this product will operate or perform substantially in conformance with the Company's published specifications and be free from defects in material and workmanship, when subjected to normal, proper and intended storage and usage by properly trained personnel, for the duration of the manufacturer's warranty period set forth in this product documentation, published specifications or package inserts. If a manufacturer's warranty period is not specified in the product documentation, published specifications or package inserts, the warranty period shall be one (1) year from the date of Bill of Sale to Buyer for equipment, and the earlier of any stated expiration dates or ninety (90) days for all other products (the "Warranty Period"). The Company agrees during the Warranty Period, to repair or replace, at the Company's option, defective Products so as to cause the same to operate in substantial conformance with said published specifications; provided that Buyer shall (a) promptly notify the

# **User Manual**

Company in writing upon the discovery of any defect, which notice shall include the product model and serial number (if applicable) and details of the warranty claim; and (b) after the Company's review, the Company will provide Buyer with service data and/or a Return Material Authorization ("RMA"), which may include biohazard decontamination procedures and other product-specific handling instructions, then, if applicable, Buyer may return the defective Products to the Company with all costs prepaid by Buyer. Replacement parts may be new or refurbished, at the election of the Company. All replaced parts shall become the property of the Company. Shipment to Buyer of repaired or replacement Products shall be made in accordance with the Delivery provisions of the Terms & Conditions agreed to by the purchase of this product. Consumables are expressly excluded from this warranty.

Notwithstanding the foregoing, Products supplied by the Company but not manufactured by the Company are not warranted by the Company, and the Company disclaims all warranties associated with the accuracy of any original manufacturer or third-party supplier Product information, but the Company agrees to assign to Buyer any warranty rights in such Product that the Company may have from the original manufacturer or third-party supplier, to the extent such assignment is allowed by such original manufacturer or third-party supplier.

In no event shall the Company have any obligation to make repairs, replacements or corrections required, in whole or in part, as the result of (i) normal wear and tear, (ii) accident, disaster or event of force majeure, (iii) misuse, fault or negligence of or by Buyer, (iv) use of the Product in a manner for which they were not designed, (v) causes external to the Product such as, but not limited to, power failure or electrical power surges, (vi) improper storage and handling of the Product or (vii) use of the Product in combination with equipment or software not supplied by the Company. If the Company determines that Product for which Buyer has requested warranty services are not covered by the warranty hereunder, Buyer shall pay or reimburse the Company for all costs of investigating and responding to such request at the Company's then prevailing time and materials rates. If the Company provides repair services or replacement parts that are not covered by this warranty, Buyer shall pay the Company therefor at the Company's then prevailing time and materials rates.

ANY INSTALLATION, MAINTENANCE, REPAIR, SERVICE, RELOCATION OR ALTERATION TO OR OF, OR OTHER TAMPERING WITH, THE PRODUCT PERFORMED BY ANY PERSON OR ENTITY OTHER THAN THE COMPANY WITHOUT THE COMPANY'S PRIOR WRITTEN APPROVAL, OR ANY USE OF REPLACEMENT PARTS NOT SUPPLIED BY THE COMPANY, SHALL IMMEDIATELY VOID AND CANCEL ALL WARRANTIES WITH RESPECT TO THE PRODUCT.

# **User Manual**

THE OBLIGATIONS CREATED BY THIS WARRANTY STATEMENT TO REPAIR OR REPLACE A DEFECTIVE PRODUCT SHALL BE THE SOLE REMEDY OF BUYER IN THE EVENT OF A DEFECTIVE PRODUCT. EXCEPT AS EXPRESSLY PROVIDED IN THIS WARRANTY STATEMENT, THE COMPANY DISCLAIMS ALL OTHER WARRANTIES, WHETHER EXPRESS OR IMPLIED, ORAL OR WRITTEN, WITH RESPECT TO THE PRODUCTS, INCLUDING WITHOUT LIMITATION ALL IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR ANY PARTICULAR PURPOSE OR USE. THE COMPANY DOES NOT WARRANT THAT THE PRODUCTS ARE ERROR-FREE OR WILL ACCOMPLISH ANY PARTICULAR RESULT.

NOTHING IN THIS MANUAL OR IN THE TERMS AND AGREEMENT AGREED TO AT THE TIME OF PURCHASE SHALL BE CONSTRUED TO PROVIDE FOR THE TRANSFER, ASSIGNMENT OR ALIENATION OF THE LIMITED WARRANTIES PROVIDED BY THE COMPANY HEREIN TO ANY OTHER PARTY OTHER THAN BUYER, INCLUDING BUT NOT LIMITED TO, SUCCESSORS, HEIRS, SUBSIDIARIES OR ANY OTHER THIRD PERSON.

THE WARRANTIES IDENTIFIED ABOVE ARE THE SOLE AND EXCLUSIVE WARRANTIES WITH RESPECT TO THIS PRODUCT AND ARE IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED. ALL OTHER WARRANTIES ARE EXPRESSLY DISCLAIMED, LIMITATION, WITHOUT THOSE OF NON-INFRINGEMENT, MERCHANTABILITY, OR FITNESS FOR A SPECIFIC PURPOSE OR USE (WHETHER ARISING FROM STATUTE, OR OTHERWISE IN LAW, OR FROM A COURSE OF DEALING, OR USAGE OF TRADE). THIS WARRANTY IS NOT TRANSFERABLE FROM THE ORIGINAL PURCHASER TO A SUBSEQUENT OWNER. FURTHER, THE COMPANY IS NOT LIABLE IN CASES OF DELIBERATE, NEGLIGENT OR ACCIDENTAL MISUSE OF THE PRODUCTS, USE WITH INAPPROPRIATE REAGENTS OR CONSUMABLES, DAMAGE CAUSED BY DISASTER, REPAIR OR MODIFICATIONS DONE BY ANYONE OTHER THAN THE COMPANY. THE COMPANY DOES NOT WARRANT THAT THE PRODUCT WILL NOT INFRINGE THE INTELLECTUAL PROPERTY RIGHTS OF A THIRD PARTY EITHER ALONE OR IN COMBINATION WITH OTHER PRODUCTS OR IN THE OPERATION OF ANY PROCESS. THE COMPANY'S TOTAL LIABILITY FOR BREACH OF THESE WARRANTIES SHALL BE LIMITED TO THE GROSS PURCHASE PRICE FOR THE PRODUCT AS STATED IN THE RELEVANT INVOICE OR ORDER UNDER WHICH THE ACTIVITY GIVING RISE TO LIABILITY ARISES. UNDER NO CIRCUMSTANCES SHALL THE COMPANY BE LIABLE FOR INDIRECT, INCIDENTAL, CONSEQUENTIAL, SPECIAL OR PUNITIVE DAMAGES (INCLUDING LOSS OF BUSINESS OR PROFITS) UNDER ANY THEORY OF LIABILITY HEREUNDER.

## **User Manual**

For claims under the warranty please contact your local supplier. You may also send the instrument directly to manufacturer, enclosing the invoice copy and by giving reasons for the claim.

#### LIMITATION OF LIABILITY

NOTWITHSTANDING ANYTHING TO THE CONTRARY CONTAINED HEREIN, THE LIABILITY OF THE COMPANY UNDER THESE TERMS AND CONDITIONS (WHETHER BY REASON OF BREACH OF CONTRACT, TORT, INDEMNIFICATION, OR OTHERWISE, BUT EXCLUDING LIABILITY OF THE COMPANY FOR BREACH OF WARRANTY (THE SOLE REMEDY FOR WHICH SHALL BE AS PROVIDED UNDER THIS SECTION) SHALL NOT EXCEED AN AMOUNT EQUAL TO THE LESSER OF (A) THE TOTAL PURCHASE PRICE THERETOFORE PAID BY BUYER TO THE COMPANY WITH RESPECT TO THE PRODUCT(S) GIVING RISE TO SUCH LIABILITY OR (B) ONE HUNDRED THOUSAND DOLLARS (\$100,000), WHICHEVER IS LESS.

With respect to any software products incorporated in or forming a part of the Product hereunder Buyer intends and agrees that such software products are being licensed and not sold, and that the words "purchase", "sell" or similar or derivative words are understood and agreed to mean "license", and that the word "Buyer" or similar or derivative words are understood and agreed to mean "licensee". Notwithstanding anything to the contrary contained herein, The COMPANY or its licensor, as the case may be, retains all rights and interest in software products provided hereunder. See the Terms & Conditions of purchase for details regarding such license.



## **IMPORTANT Safety Items to be Observed.**

## **Safety Warnings**

#### ! DANGER (may cause serious damage to property and or casualties)

- Please carefully read this User Manual prior to operating the instrument and observe the instructions on safe operation.
- Ensure that only trained staff operate this instrument
- The power source must be grounded reliably and away from any sources of electromagnetic interference.
- Ensure that the instrument and its accessories are free of any potential defects.
- Confirm that the voltage and frequency of the power supply matches the specified voltage and frequency specified prior to use.
- Ensure the cover of the junction box is properly installed prior to operation.
- Always wear applicable PPE whilst operating the instrument.
- High Temperature Hazard; Do NOT touch the pump whilst in operation
- Avoid inhalation and contact with toxic liquids, gases, spatters, vapor, dust, biological or microbiological mediums as these may be harmful to the operator.
- Place the instrument on a stable, dry, clean, antiskid, and fireproof surface.
- Ensure there is sufficient ventilation gaps and spacing above and around the instrument
- Ensure the instrument is stable and level and will not move due to vibration upon start-up.
- The instrument must always be monitored whilst in operation.
- Gases, vapours, or other substances at the exhaust may be harmful.
- When moving the Pump ensure the Pump is switched 'OFF' and the Power Supply is disconnected.
- Ensure the lifting eye is securely tightened prior to lifting the pump.
- Ensure the pump is kept in normal operating orientation during transport. Failure to do so can cause pump damage and or oil spillage.

#### ! WARNING (may cause property damage or personal injury)

- The working environment must be free of any flammable, explosive, volatile or corrosive substances. If long term storage is required, clean the pump and empty of all oil. Store the pump in a dry clean area at room temperature
- Ensure handling and disposal of used oil and other parts are compliant to local environmental laws and regulations
- Do not block the motor vent.
- Only accessories listed as 'optional accessories' for use, can guarantee operational safety.
- Ensure handling and disposal of packaging is compliant to local environmental laws and

# **User Manual**

regulations.

#### ! ATTENTION (may affect operational performance or service life)

- The instrument is deenergised only when the power cable is disconnected from the power source.
- Ensure the instrument and its accessories are protected external vibration.
- The Company reserves the right to modify the design and technical data of the pump without notice.
- Only accredited and qualified professional repair technicians can open the instrument or conduct required repairs. Persons performing repairs on the instrument other than those selected or approved by the Company shall operate to void any warranty contained hereinabove for the product





## 2.1 After Sales Support

If problems are encountered or technical support is required when installing or using the instrument, please contact serviceusa@hollandgreenscience.com

The company may provide technical assistance and information regarding the instrument or equipment or service without charge at its sole discretion. Buyer assumes sole responsibility for any reliance on or use of such assistance and information, and the company makes no warranty thereon.

Upon contact the following information is required:


- Product serial number (located on the instrument nameplate)
- Description of issue or problem
- Method and or operating steps you have undertaken towards resolution.
- Your contact details inclusive of telephone number and email address.





# III. Aether 8 Eco Vacuum Pump

## 3.1 Components







## 3.2 Before using the vacuum pump

- The vacuum pump's AC motor requires 110 V rms  $\pm$  10 % .
- The voltage and frequency of the power outlet should match the specifications listed on the motor's decal. The ON-OFF switch should be in the OFF position before plugging the pump into a power outlet. Remove and discard the exhaust plug from the end of the pump's handle prior to operation.
- The pump is shipped without oil in the reservoir. Before operating the pump, remove the OIL FILL cap and fill the pump with oil until the oil level reaches the bottom of the viewing glass.
   The oil capacity of the pump is approximately 340 ml. After filling the reservoir, replace the OIL FILL cap.
- Remove the cap from one of the air inlet ports and then turn the motor switch to ON. When
  the pump starts to run smoothly (may take from 2 to 30 seconds depending on the ambient
  temperature), replace the cap on the air inlet port. After the pump runs for about 1 minute,
  check the viewing glass for the oil level. The oil level should be at the viewing glass OIL
  LEVEL line. Add more oil if needed.
- Note: the oil level should remain at the OIL LEVEL line when the pump is running.
   Underfilling the pump will cause inefficient operation and reduce the level of vacuum achievable. Overfilling the pump may result in oil in the pump's exhaust.

## 3.3 Post Processing Vacuum Pump Maintenance

After using the vacuum pump, it is recommended that you do the following steps to extend the lifespan of the pump:

- Close the manifold valve between the pump and the system.
- Remove the vacuum hose from the pump inlet.
- Put a cap on the inlet port to prevent any contamination or loose particles from entering the port.




# IV. Technical Specifications

# 4.1 Specifications

Model		10611007
Product Name		Aether 8 Eco
Voltage		110V/60Hz
Free Air Displacement	CFM	5.0
	L/min	142
Ultimate Vacuum	Pa	3 x 10·1
	Mbar	0.003
	Microns	25
Moto	1/2	
Intake fit	1/4" & 3/8" Flare	
Oil Capa	340	
Dimension	340x138x244	
Net Weigl	11.4	



# 4.2 Expanded Drawing



4	0	00	T(1
1	Cross screw	20	Trestle cover board
2	Fan cover	21	Oil fill port
3	Fan	22	Seal
4	Motor cover	23	Inlet fitting
5	Bearing	24	Trestle
6	Cross screw	25	Pump 25-1 Rotary-vane
			body 25-2 Rotary-vane spring
7	Motor rotor	26	Cap board
8	Power supply cords	27	Exhaust fitting
9	Power switch	28	Die cast aluminum housing
10	Bearing	29	Viewing glass
11	Centrifugal switch	30	Oil cap
12	Insulating bushing	31	Seal
13	Handle slipcover	32	Screw
14	Nut	33	Seal
15	Handle	34	Screw
16	Capacitor	35	Rubber foot
17	Junction box	36	Screw
18	Motor hull	37	Base
19	Screw		



## VI. Proper Use

The instrument is designed for non-residential use and to be used in conjunction only with accessories recommended within this manual and by the manufacturer.


## VII. Inspection

#### Packing List

Unpack the equipment carefully and check for any damage which may have arisen during transport. In the event of identified damage, please contact serviceusa@hollandgreenscience.com

The package includes the following items

Item Description	Quantity
Main Unit	1
Quick Setup Card	1





#### VIII. Maintenance

## 8.1 Vacuum pump oil

The type of oil used in the vacuum pump determines the ultimate vacuum the pump can reach. We recommend using "High Vacuum Pump Oil" for this vacuum pump. This oil is designed to maintain maximum viscosity at normal operation temperature. This oil is also suitable for cold weather conditions.

### Oil change procedure

- Turn the vacuum pump on for several minutes and then turn it off.
- Remove the OIL DRAIN cap and drain the used oil into a container and dispose of it
  properly. Oil can be forced from the pump by opening the inlet and partially blocking the
  exhaust with a cloth while the pump is running. Do not operate the pump this way for more
  than 20 seconds.
- When the oil stops flowing out, tilt the pump forward to drain the remainder of the oil.
- Replace the OIL DRAIN cap. Remove the OIL FILL cap and fill the reservoir with new vacuum pump oil until the oil level reaches the bottom of the viewing glass. The oil capacity of the pump is approximately 340 ml.
- Check and make sure that the inlet ports are capped before turning on the pump. Let the
  pump run for one minute and then check the oil level. If the oil level is below the OIL LEVEL
  line on the viewing glass, add more oil slowly (with the pump running) until the oil level
  reaches the OIL LEVEL line. Replace the OIL FILL cap. Make sure that the inlet is capped,
  and the drain cap is on tight.
- If the oil is contaminated with sludge that forms when water is allowed to collect in the oil, remove the oil reservoir cover and wipe it out.
- Another way of dealing with contaminated oil is to force out the oil from the pump reservoir. To do this, let the pump run until it is warmed up. While the pump is still running, remove the oil drain cap and slightly restrict the exhaust. This will back-pressure the oil reservoir and force out the contaminated oil. When the oil stops flowing out, turn off the pump. Repeat this procedure until the contaminated oil is removed. Replace the OIL DRAIN cap and refill the reservoir with fresh pump oil to the proper level.

## **User Manual**

## IX. Troubleshooting

If dissembling of the pump is required, please check your warranty first to ensure that the warranty will not be voided due to dissembling by unauthorized service personnel.

#### 9.1 Failure to start

Check the line voltage. The vacuum pump requires 110 Vrms ±10% (loaded) at 32 °F.
 Operating outside the recommended range may result in damage to the system.

## 9.2 Oil leakage

 The module cover gasket or the shaft seal may need to be replaced. If leakage occurs in the oil drain plug, reseal the plug using a commercial pipe thread sealer.

### 9.3 Failure to attain vacuum

- Make sure the connection of vacuum gauge in good. Leakage can be confirmed by
  monitoring the vacuum with a thermistor gauge while applying vacuum pump oil at
  suspected leak points. The vacuum will improve briefly when the oil is sealing the leak.
- Make sure the pump oil is clean. A contaminated pump may require several oil flushes.
- For optimal pump efficiency, the oil level should reach the OIL LEVEL line on the viewing
  glass when the pump is running. Do not overfill. The heat generated while the pump is
  running will cause the oil to expand and appear at a higher level than when the pump is not
  running. To check the oil level, start the pump with the inlet capped. Check the oil level
  through the viewing glass. Add more oil if needed.