

XIROS MIKRO V1.12 FREEZE DRYER: COMPLETE PRESERVATION BUNDLES



WE ARE EXCITED TO OFFER OUR NEW PRESERVATION BUNDLES, DESIGNED FOR THOSE WHO ARE FOCUSED ON FREEZE DRYING CANDY, FRUITS AND OTHER CONSUMABLE SNACKS. EACH BUNDLE INCLUDES THE SKUS LISTED BELOW FROM HOLLAND GREEN SCIENCE

ONE XIROS MIKRO V1.12 FREEZE DRYER



THE XIROS MIKRO SMALL FREEZE DRYER V1.12 IS SOLD STANDARD AS A NINE-SHELF FREEZE DRYER IN THE UNITED STATES. HOWEVER, THE XIROS MIKRO HAS THE ABILITY TO EASILY BE CONVERTED TO A THREE, FIVE OR SEVEN SHELF ASSEMBLY THANKS TO OUR INTERCHANGEABLE SHELVING SYSTEMS.


ONE ANEMOS DRY SCROLL PUMP

THE ANEMOS DRY SCROLL VACUUM PUMP IS THE LATEST DESIGN IN OIL-FREE VACUUM PUMP TECHNOLOGY. THIS PUMP PROVIDES OIL AND HYDROCARBON FREE VACUUM TO PREVENT CONTAMINATION. THE ANEMOS IS A GREAT REPLACEMENT FOR ALL ROTARY VANE PUMPS AND CAN BE USED IN MULTIPLE APPLICATIONS INCLUDING FREEZE DRYING. NO MAINTENANCE REQUIRED AND DESIGNED TO RUN 24/7.



ONE 3, 5 OR 7 SHELF STACK + TRAYS

SHELF STACKS ARE AVAILABLE IN 3, 5 OR 7 SHELF ASSEMBLY OPTIONS. SIMPLY DISCONNECT THE 9 SHELF SYSTEM FROM YOUR XIROS MIKRO, REMOVE AND REPLACE WITH THE SHELF STACK OF YOUR CHOOSING. TRAYS ARE ALSO INCLUDED FOR EACH ASSEMBLY.



7 SHELF BUNDLE

SHELF SPACE = 22mm
TRAY HEIGHT = 15mm
PRICE = \$10,650 USD

5 SHELF BUNDLE

SHELF SPACE = 33mm
TRAY HEIGHT = 20mm
PRICE = \$10,500 USD

3 SHELF BUNDLE

SHELF SPACE = 59mm
TRAY HEIGHT = 20mm
PRICE = \$10,350 USD

CONTACT US

SALES@HOLLANDGREENSCIENCE.COM

US.HOLLANDGREENSCIENCE.COM

HOLLAND GREEN SCIENCE 1440 W. TAYLOR ST # 1817 CHICAGO, IL 60607